



# Multiplex detection of KRAS mutations in colorectal cancer FFPE samples using droplet digital PCR



Wei Yang, Dawne N Shelton, Jennifer R Berman, Bin Zhang, Samantha Cooper, Svilen Tzonev, Eli Hefner, John F Regan  
Digital Biology Center, Bio-Rad Laboratories, 5731 W. Las Positas Blvd, Pleasanton, CA, 94588

### Introduction

Targeted therapies in many cancers have allowed unprecedented progress in the treatment of disease. However, routine implementation of genomic testing is limited due to: 1) difficulties in detection of mutational loads below 5%, 2) limited amounts of sample (pg-ng range) per biological specimen, 3) diagnostic turnaround time and workflow, and 4) cost. KRAS is mutated in approximately 40% of colorectal cancers, and KRAS mutations are indicative of a negative response to cetuximab therapy. To optimize therapy strategies for personalized care, it is therefore critical to rapidly screen patient samples for the presence of multiple KRAS mutations. We have developed a multiplexing strategy to screen clinically-actionable KRAS mutations using digital PCR. No pre-amplification step was required. This sensitive and inexpensive method reduces the risk of contamination and can be easily implemented in molecular diagnostic laboratories for rapid, routine screening of cancer patients.

### Materials and Methods



- 16 mCRC (7 female, 9 male, average age 64 years) and 4 grossly normal colon (2 female, 2 male, average age 65 years) FFPE blocks were purchased (Advanced Tissue Services, Phoenix, AZ). mCRC samples were classified as KRAS-mutation positive by the vendor. Samples were prepared using standard protocols (Qiagen).
- ddPCR was performed on 1-5 µl per sample per well using either a multiplexed KRAS G12/G13D assay or validated PrimqPCR ddPCR mutation assays for one of seven individual KRAS mutations (G12D, G12V, G13D, G12A, G12C, G12R, G12S, Bio-Rad Laboratories)
- Positive mutation references were from Horizon Diagnostics, and negative controls were wildtype-only from Promega (Female gDNA). Statistical significance was determined using 95% confidence intervals.

Figure 1: Multiplexed single-well detection of 7 actionable KRAS mutations




Figure 2: FFPE samples yield low and variable amounts of amplifiable DNA




Figure 3: Multiplex detection of KRAS mutations in 16 mCRC FFPE samples




Figure 4: Correlation of quantified KRAS mutational load by multiplex and single KRAS assays




Figure 5: Sample heterogeneity is easily detected




Figure 6: KRAS mutations detected in grossly normal FFPE samples




Figure 7: Multiplex screening of additional actionable EGFR pathway mutations in a single well



### Conclusions

- We have demonstrated sensitive and precise detection (less than 1%, single reaction) of multiple actionable KRAS mutations in FFPE samples from colorectal cancer patients
- Concordance between duplex- and multiplex-based detection is excellent
- Precision and sensitivity of the measurements reveal significant spatial heterogeneity of the samples
- Droplet-digital PCR provides a simple and robust workflow for mutation detection in clinical patient samples
- Future work will extend the number of targets per reaction, optimize the workflow and sample compatibility to deliver on the promise of NextGen Digital technologies in precision medicine and research