

# Quantification of 100 RNA targets from limited samples or single cells

Steven Okino, Michelle Kong, and Yan Wang

Gene Expression Division, Life Science Group, Bio-Rad Laboratories, Inc., 2000 Alfred Nobel Drive, Hercules, California, 94547, USA



## Gene expression profiling during differentiation

# Single cell analysis

#### Abstract

Gene expression profiling is often limited by small sample size and low target gene expression. To address these issues we developed SsoAdvanced<sup>™</sup> PreAmp Supermix, a PCR-based reagent that increases the expression of up to 100 target genes at least 1000-fold. To validate our reagent, we analyzed a panel of genes that are regulated during stem cell differentiation. We find that relative gene expression results generated using this product are statistically equivalent to results obtained using standard qPCR; however, 10,000-fold less sample can be used. Importantly, SsoAdvanced PreAmp Supermix maintains patterns of gene expression changes across samples, thus the same biological insights would be derived from a pre-amplification (PreAmp) experiment and a standard gene expression profiling experiment. In addition, we demonstrate the application of this product for analyzing single cell gene expression.

#### Table of genes analyzed

We analyzed a panel of 100 genes that are involved in stem cell differentiation. PrimePCR<sup>™</sup> PreAmp and qPCR assays were used in the PreAmp reaction and for qPCR analysis, respectively.

Panel of 100 Genes Analyzed											
gDNA	POU5F1	ZFP42	CD44	DPPA5	GATA2	KLF4	NK X2-5	SMAD1	TAL1		
HPRT1	POU3F2	DNMT3B	CDC42	EN2	GATA6	LIN28A	NROB1	SMAD2	TAT		
B2M	NEUROD1	NR5A2	CDK1	ENG	GDF3	MEIS1	NT5E	SMAD3	TCF3		
GAPDH	LEFTY2	NR6A1	CHD1	ESRRB	GFAP	MESP1	OLIG2	SOX15	TCL1A		
TFRC	UTF1	SOX2	CHD7	ETV2	GSC	MIXL1	OTX2	SOX17	TEK		
ТВР	TDGF1	ACTA2	CNOT3	FGF2	HAND1	MYBL2	PAF1	SOX3	TERT		
NANOG	NEUROG2	AICDA	DES	FGF5	HNF4A	MYC	RIF1	SOX7	THAP11		
PAX6	ASCL1	ALB	DPPA2	FLII	HSPA9	MYOD1	ALAS1	STAT3	THY1		
TB X3	ASCL2	ALPL	DPPA3	FOXG1	ISL1	NCAM1	SALL4	Т	TRIM 28		
PRDM14	KLF2	CCNA2	DPPA4	GATA1	KAT5	NES	SCN1A	TAGLN	ZFX		

#### Analysis of cell populations undergoing differentiation

NTera2 cells (NT2) are a well-established model system used to study human stem cell behavior<sup>1</sup>. When treated with retinoic acid (RA), NT2 differentiate into neurons<sup>2</sup>. As NT2 differentiate, the level of stem cellspecific biomarkers decrease and the level of neuron-specific biomarkers increase<sup>3</sup>. We used the NT2 model system to determine if PreAmp can accurately quantify gene expression changes in limited samples.



#### Analysis of single cells undergoing differentiation

Most studies that address how cells change are based on population analyses; understanding of how individual cells behave in a changing population is important. For example, as a NT2 population moves from a pluripotent to a neuronal state, do individual cells change in unison, or at varying rates?



### Single cell analysis

We developed a protocol that uses SsoAdvanced PreAmp Supermix to quantify gene expression in individual NT2 cells at various time points following RA treatment. We find that cells at the 0 day time point express primarily stem cell biomarkers and cells at the 9 day time point express mainly neural cell biomarkers. Interestingly, at the 3 day time point we observe that the gene expression profiles are heterogeneous - some cells are similar to the 0 day cells while others are similar to the 9 day RA treated cells. These findings imply that NT2 cells do not undergo synchronized differentiation; instead the differentiation process is heterogeneous and cell dependent. In addition, the average gene expression profiles from the individual cells correlate with the gene expression profiles from the population analysis; this validates our single cell results.



Life Science Group 2000 Alfred Nobel Drive Hercules, CA 94547 USA

#### **Pre-amplification workflow**



#### **SsoAdvanced PreAmp Supermix generates minimal bias**

To assess bias introduced by PreAmp, we compared target gene levels in a cDNA sample before and after PreAmp. For each target, we calculated the bias introduced by PreAmp and plotted it against the no PreAmp Cq values for the target genes. For 98% of genes with a no PreAmp Cq<35, PreAmp did not introduce bias by more than 0.75 Cq.



#### **Assessment of PreAmp bias**

NT2 cells were treated with RA for 0 to 7 days; RNA was isolated and converted to cDNA. Varied amounts of cDNA (100 pg, 1 ng and 10 ng) were pre-amplified with SsoAdvanced PreAmp Supermix and all 100 target genes were subsequently analyzed by qPCR. In parallel, as a no-PreAmp control, 1 μg cDNA was used to quantify the 100 gene targets directly by qPCR (10 ng/target). The experiment was performed three times independently.

The level of expression of each target gene was determined and expressed as "RNA copies per cell" after normalizing to TBP with the assumption that there are 10 TBP copies per cell. Error bars represent standard deviation of the three independent experiments. Data for target genes that show significant differentiationinduced changes are shown.

We find that, relative to standard gene expression analysis, PreAmp does not cause statistically significant bias in the gene expression profiles of all target genes analyzed, even when 10,000-fold less cDNA is used. These findings demonstrate that SsoAdvanced PreAmp Supermix provides accurate and reproducible quantification of both stem cell- and neural cell-specific biomarkers and does not bias gene expression profiling results.



per panel per panel per panel per panel



#### Table of gene expression in differentiating single cells

<u>RNA copie</u>	es per cell								
RA		TBP	OCT4	NANOG	TDGF1	PRDM14	OTX2	MEIS1	PAX6
No cells		0	0	0	0	0	0	0	0
None	Cell 1	25	131	66	69	9	9	2	0
None	Cell 2	10	443	115	84	11	8	0	0
None	Cell 3	23	287	104	29	25	52	0	1
None	Cell 4	14	368	147	105	19	67	0	0
No cells		0	0	0	0	0	0	0	0
3d	Cell 5	9	161	19	4	0	2	18	0
3d	Cell 6	14	8	8	1	0	9	4	0
3d	Cell 7	35	67	19	7	0	6	49	1
3d	Cell 8	25	5	7	14	0	11	3	0
3d	Cell 9	25	448	184	34	7	8	6	0
3d	Cell 10	18	330	105	122	2	1	4	0
3d	Cell 11	22	43	13	7	0	2	15	2
No cells		0	0	0	0	0	0	0	0
9d	Cell 12	6	0	0	1	0	0	5	2
9d	Cell 13	18	0	1	3	0	0	4	14
9d	Cell 14	10	0	0	0	0	0	4	6
9d	Cell 15	6	3	1	2	0	0	2	4
9d	Cell 16	10	0	5	8	0	0	5	9
9d	Cell 17	19	0	5	5	0	0	9	0

#### **SsoAdvanced PreAmp Supermix has a wide dynamic range**

PreAmp was performed on 10 pg to 1 µg cDNA. Selected target genes with varied expression levels were analyzed by qPCR. The results (efficiency, R^2) show that a wide linear dynamic range of cDNA input can be achieved regardless of target expression level.



#### **SsoAdvanced PreAmp Supermix is efficient for up to 24 PreAmp cycles**

PreAmp was performed on 10 ng cDNA for 12-24 cycles. qPCR traces for TBP and KLF4 are shown. Efficiency and R<sup>2</sup> values show that this product maintains exponential amplification for up to 24 PreAmp cycles.



#### Graphical representation of gene expression in differentiating single cells



#### Summary and conclusions

We have developed a new reagent, SsoAdvanced PreAmp Supermix, which allows for accurate quantification of up to 100 target genes from very small samples including single cells. We envision that this product can benefit researchers who work with limited or rare samples, and can lead to advances in the single cell analysis field.

#### References

1. Andrews PW (1984). Retinoic acid induces neuronal differentiation of a cloned human embryonal carcinoma cell line in vitro. Dev. Biol 103(2):285-293.

2. Pleasure, S. J., and Lee, V. M. (1993). NTera 2 cells: a human cell line which displays characteristics expected of a human committed neuronal progenitor cell, J Neurosci Res 35, 585-602.

3. Deb-Rinker, P., Ly, D., Jezierski, A., Sikorska, M., and Walker, P. R. (2005). Sequential DNA





