



Simultaneous Quantification of HDR and NHEJ Editing Events Induced by Site-Specific Nucleases Using ddPCR

BIO-RAD

Jennifer R. Berman¹, Samantha B. Cooper¹, Bin Zhang¹, Yuichiro Miyaoka², Amanda H. Chan², Steven J. Mayer^{1,2}, George Karlin-Neumann¹, Bruce R. Conklin²

¹Digital Biology Center, Bio-Rad Laboratories, 5731 W. Las Positas Blvd., Pleasanton, CA, 94588, USA


²Gladstone Institute of Cardiovascular Disease, 1650 Owens St., San Francisco, CA, 94158, USA

Introduction

Genome editing tools such as TALENs and the CRISPR/Cas9 system have revolutionized our ability to edit the genome of any cell including human induced pluripotent stem cells (iPSCs). Sequence-specific nucleases induce double strand breaks or nicks at target sites, activating the DNA repair pathways of non-homologous end joining (NHEJ) or homology-directed repair (HDR). Through the production of small insertions or deletions, error-prone NHEJ is useful for disrupting gene function. However, for many applications, HDR is more desirable than NHEJ, since HDR utilizes homologous donor DNA to produce precise gene repair, while NHEJ causes unpredictable damage. Unfortunately, cell populations edited by HDR often contain alleles damaged by NHEJ. Since NHEJ and HDR involve different repair enzymes, it is conceivable that conditions could be achieved with high HDR and low NHEJ. However, balance between NHEJ and HDR is elusive, since we lack a rapid sensitive assay to quantify NHEJ and HDR at endogenous genomic loci. To overcome this hurdle, we developed a method to detect NHEJ and HDR simultaneously based on droplet digital PCR (ddPCR). ddPCR is a digital PCR technology based on water-oil emulsions to partition a reaction into over 20,000 uniform droplets. Combination of ddPCR and fluorescent oligonucleotide probes specific to wild-type, NHEJ, and HDR alleles allows us to simultaneously detect NHEJ and HDR events induced by site-specific nucleases in a sensitive and quantitative manner. By using our ddPCR-based method, we are able to test multiple conditions for genome editing including different types of sequence-specific nucleases and donor DNAs in order to tilt the balance between NHEJ and HDR towards HDR. We are currently optimizing mutagenesis conditions to activate HDR with minimum NHEJ in human iPSCs. Our system will result in significant improvements to genome editing technology.


Materials and Methods

QX200 Droplet Digital PCR System



- Three disease-relevant loci (RBM20, AITP78, GRN) were targeted for editing using TALEN, WT-Cas9, dCas9-FokI, or Cas9-Nickase strategies. TALENs were constructed using the Voyta laboratory's Golden Gate assembly system (Cormak et al., *Nucleic Acids Res.* 39, e82, 2011) and the MDO15 backbone vector (gift from M. Porteus and M. Rahdar, Stanford University). Guide RNAs (gRNAs) were cloned into the Zhang laboratory's pX335 (nickase) and pX330 (nuclease) (Cong et al., *Science*, 819-23, 2013; Ran et al., *Cell* 154, 1380-1389, 2013).
- QX100 & QX200 Droplet Digital PCR systems were used (Bio-Rad). ddPCR reactions were assembled using standard protocols: ddPCR Supermix for Probes (no dUTP) was combined with 5 ng – 125 ng sample, 1 μ L of 20x assay (each 1x primer at 900 nM, and if necessary a dark non-extendible probe at 500 nM, all from IDT), 4 units of restriction enzyme (HaeIII for AITP78 and RBM20 or MseI for GRN targets), and water to 20 μ L. Reactions were converted into ~20,000 1-nL droplets using the Droplet Generator and transferred to a 96-well plate for thermal cycling. After thermal cycling, droplets were read on the Droplet Reader and assigned as positive or negative based on fluorescence amplitude.

A. Targeting design for RBM20 & GRN loci



- (A) Targeting strategy for introducing pathological point mutations in human induced pluripotent stem (iPS) cells and HEK293 cells for disease-relevant loci. A pair of gRNAs (shown in magenta) for Cas9 nickase, CRISPR-Cas9 (singly) and dCas9-FokI as well as TALENs were designed to introduce the disease-relevant point mutation in the RBM20, AITP78 (not shown), and GRN loci. The wild-type and mutant residues are highlighted by a blue and a red rectangles, respectively. The nick/cut sites are indicated by red triangles.
- (B) Sequences of an isolated GRN R989X iPS cell clone aligned to the WT sequence. This clone had the C>T point mutation and a 24-bp deletion. These results indicate that HDR and NHEJ can occur in single cells, resulting in compound heterozygous mutant cells.

Figure 1: ddPCR assay for simultaneous single-well quantification of HDR and NHEJ





Figure 2: Design and Validation of RBM20 R636S Total Edit Assays



- Quantify the relative contributions of HDR and NHEJ alleles under different editing conditions, 4 Total Edit assays were designed for RBM20 R636S editing: CRISPR-Cas9, Cas9-Nickase, TALEN, and dCas9-FokI (not shown). A schematic of probe position is shown for each assay, and a representative 2D plot on a WT-only and positive control (gBlock) sample.
- Limit of Detection (LoD) for HDR and NHEJ alleles was empirically determined by assaying 100ng WT DNA (Promega) plus a 2-fold serial dilution of globs containing the point mutation edit (HDR positive control) or a 1-bp deletion at the cut site (NHEJ positive control). LoD was <0.05% for HDR and <0.1% for NHEJ. Data merged from 2 technical replicate wells is shown. 95% confidence interval is displayed. The dotted lines represent the top of the 95% CI for the WT-only negative control well.

Figure 3: Absolute quantification of HDR and NHEJ alleles in edited samples




Figure 4: Differential NHEJ production at paired gRNA sites in Nickase-edited samples





Figure 5: Additional ddPCR assay strategies for counting HDR edits & NHEJ alleles



Conclusions

- Droplet digital PCR is a high-throughput, ultra-sensitive method for rapid and inexpensive quantification of genome editing events.
- Using the Total Edit Assay, HDR (targeted edit), NHEJ (indel), and WT alleles can be counted simultaneously in a single well. NHEJ and HDR can also be counted independently with a Drop-off (NHEJ) or Edit Detection (HDR) assay. The NHEJ assay provides a quantitative readout of nuclease activity, and could provide a readout of gRNA efficacy.
- The ddPCR Total Edit Assay offers a rapid readout for identification of cell pools with a high HDR:NHEJ ratio, and for optimization of genome editing protocols.

Contact: jen_berman@bio-rad.com

