



TPN
BE ON TARGET

Defining the New Normal in Quantitative Western Blot Data

By Geetha Yadav and Kenneth Oh

Do you know that journals like the JBC strongly caution against the use of housekeeping proteins for normalization? Most journals now demand rigorous standards for presentation of western blotting data. While housekeeping proteins (HKP) have traditionally been used for decades, the problems with their usage have come to the fore only recently. More validation is required for presenting data using HKP. Normalization using total proteins present in the sample has gained significant traction and has proved to be a reliable method for normalization of western blots. Here we illustrate why the **total protein normalization (TPN)** method is superior and the way to go for western blotting.

BIO-RAD

DEFINING THE NEW NORMAL

IN QUANTITATIVE WESTERN BLOT DATA

HOW DO I ACCURATELY QUANTITATE THE CHANGES IN MY PROTEINS?

OPTION
A

HKP

NORMALIZE TO HOUSEKEEPING PROTEIN (HKP)



Normalizing with HKP: like trying to hit a moving target.

The expression level of HKPs can change due to:


- Experimental conditions
- Developmental changes
- Post-transcriptional regulation
- Differences in cell, tissue age, and type
- Stages in cell cycle

OPTION
B

TPN

USE TOTAL PROTEIN NORMALIZATION (TPN) 

 **PREFERRED METHOD**



Normalizing with total protein: multiple arrows covering the entire target area.

TPN takes into account:

- The intensity of all proteins in the lane
- Sample loading variations
- Variations during electrophoresis
- Variations during transfer

TPN PROVIDES ACCURATE QUANTITATION OF TARGET PROTEINS.



Total protein normalization is the most reliable method of reporting quantitative western blot data.



TPN is the preferred method of normalization for:




- Stem cell research
- Neuroscience research
- Alzheimer's research
- Developmental biology
- Aging studies



SEVERAL GROUPS HAVE EXPRESSED CONCERNS REGARDING WESTERN BLOT DATA NORMALIZATION:



- Nature
- Science
- NIH workshops and meetings




JOURNALS LIKE JOURNAL OF BIOLOGICAL CHEMISTRY PREFER TPN FOR SUBMITTING WESTERN BLOT DATA FOR PUBLICATION.

JBC editors recommend TPN for normalizing western blots.

It is typically better to normalize western blots using total protein loading as the denominator.

We prefer that signal intensities are normalized to total protein by staining membranes with Coomassie Blue, Ponceau S, or other protein stains, and we strongly caution against the use of housekeeping proteins for normalization, unless there is clear demonstration that expression of the housekeeping protein is unaffected by experimental treatments.

Fosang AJ and Colbran RJ (2015). Transparency is the key to quality. J Biol Chem 290, 29,692–29,694.



<http://www.jbc.org/site/misc/fora.xhtml>

BIO-RAD'S V3 Western Workflow™ INCORPORATES TPN FOR NORMALIZATION OF WESTERN BLOT DATA.


V3 Western Workflow includes:

- Stain-free visualization of total proteins on gel after electrophoresis
- Quantitation of total proteins in a lane after transfer
- Normalization of bands to total proteins by specialized software



MORE AND MORE SCIENTISTS USE V3 Western Workflow AND TPN FOR DATA NORMALIZATION.

The number of papers using the V3 Western Workflow and TPN has been consistently increasing over the years.



TOTAL PROTEIN NORMALIZATION IS THE IDEAL WAY FOR ACCURATE WESTERN BLOT NORMALIZATION.

PREFERRED ★ METHOD



Learn more

about TPN and Bio-Rad's V3 Western Workflow.



Bioradiations.com is an online magazine that offers scientific professionals around the globe timely tips, techniques, and topics related to Bio-Rad products and services. The site offers new product information, technical and application notes describing experiments performed using Bio-Rad products, interviews with renowned researchers in the field, and product demonstrations in the form of system tours, videos, podcasts, etc. Visit us at bioradiations.com.

BIO-RAD

**Bio-Rad
Laboratories, Inc.**

Life Science
Group

Web site bio-rad.com **USA** 1 800 424 6723 **Australia** 61 2 9914 2800 **Austria** 43 1 877 89 01 177 **Belgium** 32 (0)3 710 53 00 **Brazil** 55 11 3065 7550 **Canada** 1 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 420 241 430 532 **Denmark** 45 44 52 10 00 **Finland** 358 09 804 22 00 **France** 33 01 47 95 69 65 **Germany** 49 89 31 884 0 **Hong Kong** 852 2789 3300 **Hungary** 36 1 459 6100 **India** 91 124 4029300 **Israel** 972 03 963 6050 **Italy** 39 02 216091 **Japan** 81 3 6361 7000 **Korea** 82 2 3473 4460 **Mexico** 52 555 488 7670 **The Netherlands** 31 (0)318 540 666 **New Zealand** 64 9 415 2280 **Norway** 47 23 38 41 30 **Poland** 48 22 331 99 99 **Portugal** 351 21 472 7700 **Russia** 7 495 721 14 04 **Singapore** 65 6415 3188 **South Africa** 27 (0) 861 246 723 **Spain** 34 91 590 5200 **Sweden** 46 08 555 12700 **Switzerland** 41 026 674 55 05 **Taiwan** 886 2 2578 7189 **Thailand** 66 2 651 8311 **United Arab Emirates** 971 4 8187300 **United Kingdom** 44 020 8328 2000

